

1) Preparation of the DNA Solution

← Extract DNA from samples with the recommended reagent and equipment.
Please use the proper DNA Extraction Method.

Recommended Reagent and Equipment:

- Reagent **DNeasy Blood & Tissue Kit (Qiagen: 69504)**
- PCR Machine Quick Bath (ThermoGen) - Sold under the GenePasQ® trademark.
GeneAmp 9700 (Abi), Veriti 200 (Abi), Mastercycler (Eppendorf), Rotor-Gene (Qiagen)



2) Preparation of Reaction Mixture


You should prepare enough volume of the working solution.

Need 9µL / test

9µL x number of tests (samples)

*1) Premix PCR Enzyme	5.0 µL
B) Primer solution	2.0 µL
*2) RNase-Free Water	2.0 µL
Total volume	9.0 µL

µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>




*1) Premix PCR Enzyme : Recommended

① Multiplex PCR Plus Kit (Qiagen: 206152 or Qiagen No. 206151)

② SYBR Premix ExTaq

*2) RNase-Free Water : Please prepare individually.

← Dispense the premix PCR solution to each PCR tube.



3) Addition of the Sample DNA Solution

← Add 1.0 µL of the sample DNA solution (Total 10 µL).

4) Amplify the DNA by PCR

SYBR Premix ExTaq		Multiplex PCR Plus Kit		Other Protocol ()	
① 95°C	2 min	① 95°C	5 min	① °C	min
② 95°C	15 sec	② 95°C	15 sec	② °C	sec
③ 60°C	15 sec	③ 60°C	30 sec	③ °C	sec
④ 72°C	15 sec	④ 72°C	15 sec	④ °C	sec
⑤ 4°C		⑤ 4°C		⑤ °C	
} 36 cycles		} 35 cycles		} cycles	

5) Prepare the Latex-working Solution during the PCR Step.

← C) Latex solution : D) Dilution buffer = 1.5 : 20


need 21.5µL / test x number of tests

C) Latex solution	1.5 µL
D) Dilution buffer	20.0 µL
Total volume	21.5 µL

µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>

* Keep at room temperature the C) Latex solution and D) Dilution buffer before using them.

← Mix the solution uniformly with a vortex mixer. (Latex-working solution)




6) Reaction using the C-PAS Strip

- ← Centrifuge the PCR tubes after PCR briefly.
- ← Add 21.5 µL of the Latex-working solution.
- ← Mix the solution uniformly by pipetting.
- ← Insert the C-PAS strips inside the PCR tubes.
- ← Leave them at room temperature for 10 minutes.
- * Perform the reaction at room temperature (20-30°C) and 40-80% humidity.

CAUTION : A low temperature and low humidity will cause FALSE-positive detections.

* Hold a part of the absorbing pad of the C-PAS with **dried hands** during use.



7) Read the PAS Line

Read the PAS Line with the aid of the GenePasQ® Data Sheet / Strip Reading Guide.