

Quick Guide Protocol

- Protocol of GenePasQ[®] DNA 1, GenePasQ[®] DNA 2, & GenePasQ[®] DNA 3

1) Preparation of the DNA solution

- ← Extract DNA from samples with recommended reagent and equipment.
Please use the proper DNA extraction method as shown below.
- ← Obtain SAMPLE [For Animal (Shrimp) sample, use **less than 10mg**] **HOMOGENIZE** manually.
- ← Add **40µl** of the **GenCheck** or **any Recommended** extraction reagent.
- ← Vortex.
- ← Centrifuge at **10000 rpm** for **30 seconds**.
- ← Heat block (**Thermo-Shaker**) at 98-100 °C for 10 minutes.
- ← **Stay on ice** for 1 minute.
- ← Centrifuge at **10000 rpm** for **5 minutes**.
- ← Transfer **5µl** of the **Supernatant liquid** to a **new** 1.5ml tube.
- ← Add **95µl** of the **RNase-free water**.



Recommended Equipment:

- PCR machine Quick Bath (ThermoGen) - Sold under the GenePasQ trademark.
GeneAmp 9700 (ABI), Veriti 200 (ABI), Mastercycler (Eppendorf)
Rotor-Gene (Qiagen)

2) Preparation of Reaction Mixture

Need 9µL / test

*1) Premix PCR Enzyme	5.0 µL
B) Primer solution	2.0 µL
*2) RNase-Free Water	2.0 µL
Total volume	9.0 µL

You should prepare enough volume of working solution
9µL x number of tests (samples)

µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>


*1) Premix PCR Enzyme : Recommended

① **Multiplex PCR Plus Kit (Qiagen: 206152)**

② **SYBR Premix ExTaq**

*2) RNase-Free Water : Please prepare individually.

- ← Dispense the premix PCR solution to each PCR tube.



3) Addition of the sample DNA solution

- ← Add 1.0 µL of the sample DNA solution (Total 10 µL).

4) Amplify the DNA by PCR

SYBR Premix ExTaq		Multiplex PCR Plus Kit	
① 95°C	2 min	① 95°C	5 min
② 95°C	15 sec	② 95°C	15 sec
③ 60°C	15 sec	③ 60°C	30 sec
④ 72°C	15 sec	④ 72°C	15 sec
⑤ 4°C		⑤ 4°C	
} 36 cycles		} 35 cycles	

(Other Protocol)		
①	°C	min
②	°C	sec
③	°C	sec
④	°C	sec
⑤	°C	
} cycles		

5) Prepare the Latex-working solution during the PCR step.

- ← C) Latex solution : D) Dilution buffer = **1.5 : 20**


need 21.5µL / test x number of tests

C) Latex solution	1.5 µL
D) Dilution buffer	20.0 µL
Total volume	21.5 µL

µL	<input type="checkbox"/>
µL	<input type="checkbox"/>
µL	<input type="checkbox"/>

- ← Mix the solution uniformly with vortex mixer. (Latex-working solution)

* Keep at room temperature the C) Latex solution and D) Dilution buffer before using them.



* Hold a part of the absorbing pad of the C-PAS with **dried hands** during use.

6) Reaction using the C-PAS

- ← Perform the following Steps in a STERILE room.
 - ← Centrifuge the PCR tubes after PCR briefly.
 - ← Add 21.5 μ L of the Latex-working solution.
 - ← Mix the solution uniformly by pipetting.
 - ← Insert the C-PAS in tube.
 - ← Leave them at room temperature for 10 minutes.
 - * Perform the reaction at room temperature (20-30°C) and 40-80% humidity.
- CAUTION** : A low temperature and low humidity will cause FALSE-positive detections.

7) Identification of the Shrimp Pathogen (s)

Read the PAS Line with the aid of the GenePasQ[®] Data Sheet / Strip Reading Guide.

All rights reserved. No part of this publication may be reproduced, distributed or transmitted in any form or by any means, including photocopying, recording, or other electronic or mechanical methods, without the prior written permission of the owner.

Copyright © 2018 SURE Marketing Company, Inc.

This GenePasQ[®] Kit is produced by TBA, Co., Ltd. and uses DNA Chromatography Strips: C-PAS (F6) with a Patent licensed from NGK Insulators, Ltd.

- Exclusive Distributor: PT. SURE Marketing Company

Ruko Piazza Venezia
(Palais de Europe)
Jl. Bulevar Eropa No.33
Lippo Karawaci, Tangerang 15115
Banten - INDONESIA
E-Mail: smci@sure-bio.com
URL: www.sure-bio.com

- The primers were developed by Tokyo University of Marine Science and Technology (TUMSAT) under Prof. Ikuo Hirono.